



Approved by:	Version: 20230906
Date:	Author: FW, XR <span style="float: right;">Page 1 by 1</span>

## Abbott Afinion AS100 Instrument-tag

The document describes the connection and labeling of the meter, to work with the Zafena POC-Workstation. Physical connection is between, RS232 serial to USB-A.  
Cable product code: **ZAF-523**



## Instrument-tag: LMC-AAF1

Print and stick the instrument tag on to or next to the instrument. Follow the documentation "Quick guide" to perform analysis and submit results.



## Communication settings

There are no configuration options on this meter for RS232 communication.